

PHILIPS

AVENT

Manual breast
pump

with VIA Storage System

SCF300/12



Express more. Store easily.

Express, store and feed in easy steps

The unique Philips AVENT Manual Breast Pump SCF300/12 offers you total control of the suction level and pumping rhythm. It assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

Benefits and features:

- Gently stimulates fast milk flow

Clinically proven effective

- Patented Let-down Massage Cushion

Fits all Philips AVENT breast pumps and teats

- Use adaptor to express and feed

For safe storage and transport

- Leak-proof twist-on system

Specifications

Country of origin

England

What is included

Manual breast pump: 1 pcs

Milk storage cup lid: 5 pcs

Milk storage cup adaptor: 2 pcs

Newborn nipple travel pack: 2 pcs

Spare parts: 1 pcs

Extra soft Newborn Flow nipple: 2 pcs

Pre-sterilized VIA Cup (180 ml/ 6 oz): 5 pcs

Development stages

Stages: 0 - 6 months

Material

BPA free*

