

# PHILIPS

## AVENT

Manual breast  
pump

with disposable breast pads

SCF300/15



## Express more. Quickly.

### Enhanced efficiency from a manual breast pump

A manual breast pump offers you total control of the suction level and pumping rhythm. The unique Philips Avent Manual Breast Pump assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump\*

#### **Clinically proven effective**

- Patented Let-down Massage Cushion

#### **Flexible and portable**

- Easy to assemble and use

#### **Helps you breastfeed longer**

- More milk, naturally

#### **Stimulate and express**

- Soft massage cushions trigger natural let-down, like baby

#### **Disposable Breast Pads**

- Unique four layered construction for maximum dryness

# Specifications

## Country of origin

England

## What is included

Manual breast pump: 1 pcs

Disposable breast pad: 20 pcs

Breast milk container (125 ml/4 oz): 1 pcs

Newborn teat travel pack: 1 pcs

Bottle stand/funnel cover: 1 pcs

Extra-soft newborn flow teat: 1 pcs

Sealing disc for milk storage: 1 pcs

## Development stages

Stages: 0–6 months

