

PHILIPS

AVENT

Manual breast
pump

SCF300/20



Express more. Quickly.

Breast pump inspired by nature

The unique Philips AVENT BPA Free Manual Breast Pump SCF300/20 assures you enhanced efficiency, being clinically proven to express more milk than a hospital grade double electric pump*.

Clinically proven effective

- Patented Let-down Massage Cushion

Stimulate and Express

- Gentle vacuum mimics baby's suckling for steady milk flow
- Soft massage cushions trigger natural let-down, like baby

Flexible and portable

- Easy to assemble and use

Benefits and features:

- Clinically proven

Highlights

Patented Let-down Massage Cushion



Five, soft petal massagers of the Philips AVENT breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

Soft massage cushions



Patented massage cushions of the Philips AVENT breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Easy to assemble and use



Quiet and portable.

Clinically proven



Clinically proven to express more milk than a hospital grade double electric pump*.

Specifications

Country of origin

England

What is included

- Newborn nipple travel pack: 1 pcs
- Bottle stand/funnel cover: 1 pcs
- Breast milk container (125 ml/ 4 oz): 1 pcs
- Extra soft Newborn Flow Nipple: 1 pcs
- Manual breast pump: 1 pcs
- Sealing disc for milk storage: 1 pcs

Development stages

Stages: 0 - 6 months

