

PHILIPS

AVENT

Single electronic
breast pump

SCF302/01



Designed for Comfort

Breast pump inspired by nature

The unique Philips Avent BPA Free Manual Breast Pump SCF302/01 assures you enhanced efficiency, being clinically proven to express more milk than a hospital grade double electric pump*.

Benefits and features:

- Learns and continues your pumping rhythm

Stimulate and Express

- Gentle vacuum mimics baby's suckling for steady milk flow
- Soft massage cushions trigger natural let-down, like baby

Helps you breastfeed longer

- More milk, naturally

Use three ways

- By mains power, battery power and manual

Highlights

Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Learns and continues your pumping rhythm



Electronic memory learns and at the touch of a button, continues your personal pumping rhythm.

More milk, naturally



Expression allows breast milk to be available, even when you can't be there

Mains, battery and manual



In addition to mains power, this single electronic pump can also operate using battery power and as a manual pump, for total flexibility.

Specifications

Power

Voltage: 100 - 240 V

Country of origin

England

What is included

Electronic Breast Pump: 1 pc

Bottle stand/funnel cover: 1 pcs

Battery pack: 1 pcs

Manual pump parts for when power is not available: 1

Spare parts: 1 pcs

Breast milk container (125 ml/ 4 oz): 1 pcs

Extra soft newborn flow nipple: 1 pcs

Newborn nipple travel pack: 1 pcs

Development stages

Stages: 0 - 6 months

