

PHILIPS

AVENT

Manual breast
pump

with milk storage cups

SCF310/12



Designed for Comfort

Optimum comfort breast pump with massage cushions

Being stressed or in a hurry can make it more difficult to release your milk and can also affect your milk supply. The Philips Avent breast pump SCF310/12 has been designed to increase your comfort when you are pumping.

Designed for comfort

- Patented soft 5-petal massage cushion

Clinically proven effective

- Clinically proven results*

Other great features

- Unique system for easy milk storage
- Easy manual operation

Highlights

Clinically proven results*

The gentle-draw vacuum of the Philips AVENT breast pump mimics a baby's natural suckling action and gets more milk than a hospital grade double electric pump*

Patented soft massage cushion



The unique active massage cushion is intended to help stimulate natural let-down

Unique milk storage system



Make your life easier by expressing directly into any of our wide range of Philips AVENT baby bottles and breast milk storage containers for the fridge or freezer

Specifications

Material

BPA free*

Development stages

Stages: 0 - 6 months, Pregnancy

What is included

- Manual breast pump: 1 pcs
- Extra soft Newborn Flow Nipple: 1 pcs
- Disposable breast pad: 1 pcs
- Milk storage cup: 2 pcs
- Spare parts: 2 pcs
- Newborn nipple travel pack: 1 pcs
- Milk storage cup adaptor: 1 pcs
- Milk storage cup lid: 2 pcs

Country of origin global

England

