

# PHILIPS

## AVENT

Manual breast  
pump

Breast pumps

with milk storage cups

SCF310/12



## Designed for Comfort

Unique active massage cushion for optimum comfort

Being stressed or in a hurry can make it more difficult to release your milk and can also affect your milk supply. The Philips AVENT breast pump SCF310/12 has been designed to increase your comfort when you are pumping.

### **Designed for comfort**

- Patented soft 5-petal massage cushion

### **Clinically proven effective**

- Clinically proven results\*

### **Other great features**

- Unique system for easy milk storage
- Easy manual operation

## Highlights

### Clinically proven results\*

The gentle-draw vacuum of the Philips AVENT breast pump mimics a baby's natural suckling action and gets more milk than a hospital grade double electric pump\*

### Patented soft massage cushion



The unique active massage cushion is intended to help stimulate natural let-down

### Unique milk storage system



Make your life easier by expressing directly into any of our wide range of Philips AVENT baby bottles and breast milk storage containers for the fridge or freezer

## Specifications

### Country of origin

England

### Material

BPA free\*

### Development stages

Stages: 0 - 6 months, Pregnancy

### What is included

Manual breast pump: 1 pcs

Extra soft Newborn Flow Nipple: 1 pcs

Disposable breast pad: 1 pcs

Milk storage cup: 2 pcs

Spare parts: 2 pcs

Newborn nipple travel pack: 1 pcs

Milk storage cup adaptor: 1 pcs

Lid: 2 pcs

