

PHILIPS

AVENT

Manual breast
pump

Includes 4oz bottle

SCF310/20



Designed for Comfort

Optimum comfort breast pump with massage cushions

Being stressed or in a hurry can make it more difficult to release your milk and can also affect your milk supply. Our Philips Avent breast pump SCF310/20 has been designed to increase your comfort when you are pumping.

Designed for comfort

- Patented soft 5-petal massage cushion

Clinically proven effective

- Clinically proven results*

Other great features

- Unique system for easy milk storage
- Easy manual operation

Highlights

Clinically proven results*



The gentle-draw vacuum of the Philips AVENT breast pump mimics a baby's natural suckling action and gets more milk than a hospital grade double electric pump*

Patented soft massage cushion



The unique active massage cushion is intended to help stimulate natural let-down

Unique milk storage system



Make your life easier by expressing directly into any of our wide range of Philips AVENT baby bottles and breast milk storage containers for the fridge or freezer

Specifications

Country of origin global

England

Country of origin for China

China

Material

BPA free*

Development stages

Stages: Pregnancy, 0 - 6 months

What is included

- Manual breast pump: 1 pcs
- Breast milk container (125 ml/ 4 oz): 1 pcs
- Bottle stand/funnel cover: 1 pcs
- Extra soft Newborn Flow Nipple: 1 pcs
- Sealing disc for milk storage: 1 pcs
- Spare parts: 2 pcs
- Newborn nipple travel pack: 1 pcs

