

# PHILIPS

## AVENT

Twin electronic  
breast pump

Breast pumps

Includes 2 4oz bottles

SCF314/02



## Breast pump designed for comfort

Unique active massage cushion for optimum comfort

Being stressed or in a hurry can make it more difficult to release your milk. So feeling comfortable and relaxed really helps. The Philips AVENT SCF314/02 breast pump has been designed to increase your comfort when you are pumping.

### **Designed for comfort**

- Remembers your personal rhythm
- Patented soft 5-petal massage cushion

### **Clinically proven effective**

- Gentle with clinically proven results\*

### **Other great features**

- Unique system for easy milk storage

# Highlights

## Clinically proven results



Our gentle electric breast pump is so effective that it's clinically proven to express similar amounts of milk as compared to a hospital grade electronic pump\*

## Patented soft massage cushion



The unique active massage cushion is intended to help stimulate natural let-down

## Remembers your personal rhythm



Only Philips AVENT's electric breast pumps have a unique memory feature that learns and continues your personal pumping rhythm for maximum comfort

## Unique milk storage system



Make your life easier by expressing directly into any of our wide range of Philips AVENT baby bottles and breast milk storage containers for the fridge or freezer

# Specifications

## Country of origin

England

## Material

BPA free\*

## Development stages

Stages: Pregnancy, 0 - 6 months

## Power

Voltage: 100 - 240 V

## What is included

- Bottle stand/funnel cover: 2 pcs
- Breast milk container (125 ml/ 4 oz): 2 pcs
- Electronic Breast Pump: 1
- Extra soft Newborn Flow Nipple: 2 pcs
- Microfiber travel bag: 1 pcs
- Spare parts: 3 pcs
- Newborn nipple travel pack: 2 pcs
- Disposable breast pad: 2 pcs

