

PHILIPS

AVENT

Comfort Manual breast pump

Natural

Includes 3 milk storage cups



SCF330/12



Milk for baby with comfort for mom

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Gently stimulates natural let down and milk flow

- Soft massage cushion with warm feel

Express, store and feed in one

- Includes 3 versatile milk storage cups

Easy to use, store and transport

- Compact lightweight design

Other benefits

- Intuitive assembly. Easy visual matching of parts
- Compatible with other Philips Avent feeding products
- Easy cleaning thanks to the small number of separate parts
- Easy manual operation



asimpleswitch.com

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. During pumping, it gently mimics your baby's suckling for natural let down and clinically proven effectiveness. The soft massage cushion is designed to comfortably fit the majority of moms. Should you need it, a larger size cushion is available separately (availability may vary per country).

Express, store and feed in one



Use the same cup to express, store and feed breast milk to your baby. Ideal for storage in fridge or freezer

Compact design



The breast pump has a compact design, which makes it easy to position. Its ergonomically shaped handle offers a comfortable hold and full control when expressing milk. Being small and lightweight means it is easy to store and transport, making pumping on-the-go more discreet.

Intuitive assembly



Easy visual matching of parts for intuitive assembly.

Fully compatible range



This breast pump can be used in combination with other feeding products in the Philips Avent range, including our Classic bottles and milk storage containers. Philips Avent also provides a range of breast care accessories to help you breastfeed for longer and enhance your comfort.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. All parts are dishwasher proof.

Easy manual operation

The manual breast pump is ideal if you express milk occasionally and value compactness. Easily operated with one hand.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Country of origin

England

Material

Breast pump: BPA free*

VIA milk storage cup: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Breast pump body with handle: 1 pcs

Standard size cushion: 1 pcs

Milk storage cups: 3 pcs

Milk storage cup adaptor: 1 pcs

Extra soft Newborn Flow Nipple: 1 pcs

Travel cover: 1 pcs

Breast pad sample packs: 2 (2 Day pads & 2 Night pads) pcs

Breast pump

Material: BPA free*

Design

Breast pump design: Compact design

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Development stages

Stages: Pregnancy, 0 - 6 months

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Includes VIA milk storage cups: Express, store and feed in one

