

# PHILIPS

## AVENT

### Manual Breast Pump

Comfort

Includes 125 ml (4 oz) Natural bottle



SCF330/40



## More comfort, more milk\*

### Clinically proven comfort\*

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

#### **Clinically proven less leaning forward\***

- Relaxing position thanks to its special design

#### **Soft massage cushion with massaging petals**

- Soft massage cushion gently stimulates milk flow

#### **Easy to set up, express and clean**

- Compact lightweight design
- Few separate parts and intuitive design

# Highlights

## Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping helps your milk to flow more easily.

## Soft massage cushion



Our massage cushion has a soft, velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk

flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let-down.

## Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on the go more discreet.

## Easy to set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

# Specifications

## Material

**Breast pump:** BPA-free\* (food contact parts only)

**Bottle:** BPA-free\*, Polypropylene

**Teat:** BPA-free\*, Silicone

## What is included

**Breast pump body with handle:** 1 pcs

**Standard massage cushion:** 1 pc (19.5 mm [0.77"], a larger size cushion is available separately)

**Milk container 125 ml (4 oz):** 1 pcs

**Travel cover:** 1 pcs

**Sealing disc for milk storage:** 1 pcs

## Design

**Breast pump design:** Compact design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Breast pump use:** Easy cleaning, Fully compatible range, Intuitive assembly

## Development stages

**Stage:** 0–6 months

## Functions

**No leaning forward:** Sit in a comfortable position

**Soft massage cushion:** Gentle stimulation



\* Clinically proven comfort: Tests carried out on 110 mothers showed a significant preference for Philips Avent vs. a leading competitor

\* More comfort: 73% of 73 UK breastfeeding mothers agreed that this breast pump was more comfortable than their current pump (main market brands).

\* More milk: Independent research has shown that there may be a link between stress levels and milk production. See [www.philips.com/Avent](http://www.philips.com/Avent)

\* BPA-free breast pump: Only associated with the bottle and other parts that come into contact with breast milk. Following EU regulation, 10/2011