

PHILIPS

AVENT

Manual Breast Pump

Ultra Comfort

Includes 125 ml (4 oz) Natural bottle



SCF330/70



More comfort, more milk*

Clinically proven comfort*

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Clinically proven less leaning forward*

- Relaxing position thanks to its special design

Soft massage cushion with massaging petals

- Soft massage cushion gently stimulates milk flow

Easy to set up, express and clean

- Compact lightweight design
- Few separate parts and intuitive design

Easy to combine breast and bottle feeding

- Includes our Natural bottle and teat for natural latch-on

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a soft, velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let-down.

Easy on-the-go



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on the go more discreet.

Easy to set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Specifications

Material

Bottle: Polypropylene

Teat: Silicone

What is included

Breast pump body with handle: 1 pcs
Standard massage cushion: 1 pc (19.5 mm [0.77"], a larger size cushion is available separately)

Natural bottle 125 ml/4 oz: 1 pcs

Extra-soft newborn flow teat: 1 pcs

Travel cover: 1 pcs

Sealing disc for milk storage: 1 pcs

Disposable breast pads: 12 pcs

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Development stages

Stage: 0–6 months

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

