

PHILIPS

AVENT

Single electric breast pump

Pump anytime, anywhere

Includes soft massage cushion

Natural bottle and nipple



SCF332/31



More comfort, more milk*. Anytime, anywhere

Single electric breast pump with massage cushion

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. Use our quiet pump anytime, anywhere, even with batteries. It's easy to set up, personalize, use and clean

Clinically proven less leaning forward*

- More comfortable pumping position due to unique design

1 stimulation mode and 3 expression settings

- Choose the setting that is effective & comfortable for you

Gently stimulates let down and milk flow

- Soft massage cushion with massaging petals

Easy to set up, express and clean

- Closed system, designed for hygienic expressing
- Few separate parts and intuitive design
- Includes our Natural bottle and nipple for natural latch on

Quiet pump for discreet expression, anytime

- Designed to let you express quietly

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle or container, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow comfortable for you.

Soft massage cushion



Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The iconic petal cushion is designed to mimic baby suckling, intended to gently help stimulate let down.

No milk spillage in the tubes



No milk spillage in the tubes.

Easy set up, express and clean



Few separate parts and intuitive design. Easily express milk for your baby, anytime.

Includes our Natural bottle



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Quiet pump



Designed for discreet expression.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck
Breast pump design: Compact design

Development stages

Stages: 0 - 6 months

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Settings: 1 Stimulation mode, 3 Expression settings

Soft massage cushion: Gentle stimulation

Material

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

Breast pump: BPA free* (food contact parts only)

What is included

Breast pump body: 1 pcs

Base unit incl. tubing: 1 pcs

Standard size cushion (19.5 mm): 1 pcs (a larger

size cushion is available and sold separately)

Natural bottle 4oz: 1 pcs

Travel cover: 1 pcs

Sealing disc for milk storage: 1 pcs

Disposable breast pads: 4 pcs

Breast pad sample packs: 1 pack (2pcs)

Country of origin

England



* Clinically proven comfort: In a test carried out on 110 mothers in the USA, UK, China and Russia in March 2016, mothers awarded an average of 8.6/10 score for Philips Avent performance on comfort.

* More milk: Independent research has shown that there may be a link between stress levels and milk production. see www.philips.com/AVENT

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011