


PHILIPS

AVENT

Comfort Double
electric breast pump

SCF334/13



More comfort, more milk

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

More milk in less time

- The perfect time saver for moms

Choose your most comfortable setting

- Features a gentle stimulation mode and 3 pumping settings

Gently stimulates let down and milk flow

- Soft massage cushion with warm feel

Other benefits

- Includes our Natural bottle and nipple for natural latch on
- Compact lightweight design, includes handy travel bag
- Intuitive assembly. Easy visual matching of parts
- Easy cleaning due to the small number of separate parts
- Bonus: Natural bottles, sealing discs, bottle brush and more



Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

Perfect time saver for moms

Simultaneous pumping at both breasts is proven more efficient and may even boost your ability to produce breast milk.*

Intuitive assembly



Easy visual matching of parts for intuitive assembly.

Compact design



The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. The breast pump comes with a handy travel bag.

Includes our Natural nipple



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Bonus items

Includes 2 Natural bottles (4oz), 2 sealing discs, 1 bottle and nipple brush, 12 day breast pads and 6 night breast pads.

Soft massage cushion

Our massage cushion has a soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. During pumping, it gently mimics your baby's suckling for natural let down and clinically proven effectiveness. The soft massage cushion is designed to comfortably fit the majority of moms. Should you need it, a smaller size cushion is available separately.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Country of origin

England

Material

Breast pump: Polypropylene

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Breast pump body: 2 pcs

Base unit incl. tubing: 1 pcs

Natural bottle 4oz: 2 pcs

Extra soft newborn flow nipple: 2 pcs

Sealing disc for milk storage: 2 pcs

Travel bag: 1 pcs

Bonus day breast pads: 12 pcs

Bonus night breast pads: 6 pcs

Bonus bottle brush: 1 pcs

Breast pad samples: 1 pair of Comfort breast pads and 1 pair of Night breast pads

Large size cushion: 2 pcs

Design

Breast pump design: Compact design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Development stages

Stages: 0 – 6 months

Functions

No leaning forward: Sit in a comfortable position

Settings: 3 Expression settings, 1 Stimulation mode

More milk in less time: The perfect time saver

Soft massage cushion: Gentle stimulation

