

PHILIPS

AVENT

ultra soft snuggle

Toy and ultra soft pacifier

0m+

Orthodontic & BPA-Free

1x snuggle & 1x 0-6M pacifier



SCF348/01



Baby's first friend

Includes ultra soft pacifier

Philips Avent ultra soft snuggle is a plush toy with ultra soft pacifier. It is soft and lightly weighted to help make babies feel secure. The plush toy and pacifier detach for easy cleaning. It's easy for parents and baby to find and hold.

Includes ultra soft pacifier

- Cuddly soft plush toy included with ultra soft pacifier
- Soft flexible shield for more comfort and fewer skin marks*
- 98% of babies accept the textured silicone nipple*

Easy to clean

- Detachable for easy cleaning
- The pacifier can be cleaned and sterilized separately

Designed with baby's comfort in mind

- Helps you and your baby find the pacifier
- Plush toy helps keep ultra soft pacifier in place
- Plush toy is compatible with all Philips Avent pacifiers

Soft, cute, and friendly

- Pick your favorite or collect all four adorable characters

Highlights

Includes ultra soft pacifier



Cuddly soft plush toy included with our ultra soft pacifier. Made of flexible silicone, it's soft on baby's delicate skin and designed to prevent skin marks.

Soft, flexible shield



To help protect delicate skin, our ultra soft flexible shield follows the shape of baby's cheeks. It leaves fewer marks and less irritation so that babies can enjoy a more comfortable soothing experience.

98% nipple acceptance



Babies know what they like! We asked moms how their little ones respond to our textured silicone nipples and 98% said that their babies accept the Philips Avent ultra soft pacifiers.

Detaches for easy cleaning




Plush toy easily detaches from ultra soft pacifier and both are easy to clean.

Easily cleaned and sterilized



Plush toy can be machine or hand washed. ultra soft pacifier can be cleaned in dishwasher, a sterilizer or in boiling water. Feel good knowing it is always clean and safe for your little one.

Easy to find



No more hunting and searching for pacifiers! The plush toy makes the pacifier easy to find.

Helps keep pacifier in place



Plush toy is soft and lightly weighted, which helps keep ultra soft pacifier close to baby. This helps keep it and the pacifier in place while keeping baby company.

Fits Philips Avent pacifiers



Plush toy is compatible with all Philips Avent Soothies and pacifiers, so you can mix and match and create the product that is right for baby.

A little buddy for baby



Plush toy is a friendly first companion for baby. Choose from four adorable characters: giraffe, monkey, elephant or seal.

Specifications

What is included

ultra soft snuggle: 1 pcs
0-6M ultra soft pacifier: 1 pcs

Safety

BPA-free

Hygiene

Easy to clean
Pacifier can be sterilized
Pacifier is dishwasher safe
Plush toy is machine washable



- * * In a 2016 USA consumer test with 112 moms, 96% said that the shield causes fewer skin marks and less skin irritation.
- * 2016-2017 US consumer tests show an average of 98% nipple acceptance of the textured Philips Avent nipple used in our 0-6m and 6-18m ultra air and ultra soft pacifiers.
- * No 1 global pacifier brand
- * For hygiene reasons, replace pacifiers after 4 weeks of use
- * Our range supports moms and babies at every development stage