

PHILIPS

AVENT

Anti-colic baby
bottle

3 Bottles

9oz/260ml

Slow flow nipple

1m+



SCF407/33



Clinically proven to reduce colic and discomfort* Designed to reduce gas

Our anti-colic bottle's Airflex venting system and textured nipple is designed to minimize feeding interruptions and discomfort. With its integrated anti-colic valve, air is vented into the bottle and away from the baby's tummy.

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- 60% less fussing at night*
- Ribbed texture prevents collapse for uninterrupted feeding
- Nipple shape designed for secure latch

Easy to use, clean and assemble


- Leak-free design
- Easy to assemble with only fewer parts
- Easy to hold
- Wide neck bottle with rounded corner for easy cleaning

Other benefits

- Compatible range from breastfeeding to cup
- This bottle is BPA free
- Full compatibility with Classic+ bottle range
- Different nipple flow rates available


Highlights

Less colic



Our clinically proven Airflex venting system is designed to reduce gas and colic. As your baby feeds, the valve integrated into the nipple flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce gas, spit-up and burping.

Less fussing



Philips Avent anti-colic bottle reduces fussing. Babies fed with Philips Avent anti-colic bottles experienced 60% less fussing at night, than babies fed with a leading competitor's vented bottle. *

Uninterrupted feeding



The nipple shape allows secure latch and the ribbed texture helps prevent collapse for uninterrupted, comfortable feeding.

Leak-free

Our anti-colic bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

Secure latch

Nipple shape designed for secure latch

Easy to assemble



Our Anti-colic bottle has fewer parts for quick and simple assembly.

Growing up with your baby



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Easy to hold



The unique bottle shape makes this bottle easy to hold and grip in any direction.

Easy to clean



With its wide neck and rounded corners, our bottle is easy to clean quickly and thoroughly. Ease of cleaning helps to ensure more hygienic feeding.

BPA free




Philips Avent anti-colic bottle is made of BPA free material (PP).

Compatible with Classic+ range

Full compatibility with Classic+ bottle range

Fits your baby's growing needs



The Philips Avent Anti-colic bottle range offers different nipple flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All nipples are available in twin packs: Newborn, Slow, Medium, Fast and Variable Flow.

Specifications

Country of origin

United Kingdom

Material

Bottle: Polypropylene, BPA free

Nipple: Silicone, BPA free

What is included

Baby bottle: 3 pcs

Bottle

Material: BPA free

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Dishwasher & microwave safe, Easy to assemble, Easy to clean, Easy to hold

Development stages

Stages: 0-12 months

Functions

Anti-colic valve: Airflex venting system designed to reduce air ingestion

Nipple: Easy latch on, Ribbed texture prevents nipple collapse, Proven anti-colic system

Ease of use: easy to clean and assemble, leak-free, 4 pieces for easy assembly

Compatibility: Full Classic+ compatibility



- * At 2 weeks of age, babies fed with a Philips Avent bottle showed less colic, and significantly less fussing at night compared to babies fed with another competitor bottle.
- * Nipple design proven to prevent nipple collapse and associated air ingestion and feeding interruptions.
- * What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying, fussing, gassiness and spit-up.