



Anti-colic nipple

- 2 pieces
- Newborn flow nipple
- 0m+

SCF421/27

Clinically proven to reduce colic and discomfort*

Designed to reduce gas

Our anti-colic nipple is part of the Airflex venting system, air is vented into the bottle and away from the baby's tummy. The ribbed texture prevents nipple collapse and reduces feeding interruptions and discomfort.

Benefits

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- 60% less fussing at night*

Uninterrupted feeding

- Nipple shape and ribbed texture prevents collapse
- Nipple shape designed for secure latch

Other benefits

- Compatible range from breastfeeding to cup
- Different nipple flow rates available
- This nipple is BPA free
- Full compatibility with Classic+ bottle range

PHILIPS


Features

Less colic



Our clinically proven Airflex venting system is designed to reduce gas and colic. As your baby feeds, the valve integrated into the nipple flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce gas, spit-up and burping.

Less fussing



Philips Avent anti-colic bottle reduces fussing. Babies fed with Philips Avent anti-colic bottles experienced 60% less fussing at night, than babies fed with a leading competitor's vented bottle. *

Uninterrupted feeding



The nipple shape allows secure latch and the ribbed texture helps prevent collapse for uninterrupted, comfortable feeding.

Secure latch

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Country of origin

United Kingdom Yes

Material

Nipple Silicone
BPA free

What is included

Nipple 2 pcs

Development stages

Stages 0-6 months

Growing up with your baby



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Fits your baby's growing needs



The Philips Avent Anti-colic bottle range offers different nipple flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All nipples are available in twin packs: Newborn, Slow, Medium, Fast and Variable Flow.

BPA free



Philips Avent anti-colic nipple is made of BPA free material (silicone).

Compatible with Classic+ range

Functions

Anti-colic valve	Helps to reduce air ingestion (it fits Airflex venting system)
Nipple	Easy latch on, Ribbed texture prevents nipple collapse, Proven anti-colic system
Ease of use	easy to clean and assemble leak-free
Compatibility	Full Classic+ compatibility

Nipple

Flow speed	Newborn flow
------------	--------------

- * At 2 weeks of age, babies fed with a Philips Avent bottle showed less colic, and significantly less fussing at night compared to babies fed with another competitor bottle.
- * Nipple design proven to prevent nipple collapse and associated air ingestion and feeding interruptions.
- * What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying, fussing, gassiness and spit-up.



data subject to change
2019, February 24

Version: 4.1.1
UPC: 00075020062222

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com