

# PHILIPS

## AVENT

### Classic+ baby bottle

4 Bottles

4oz/125ml

Newborn flow nipple

0m+



SCF560/47



## Clinically proven to reduce colic and discomfort\*

### Few parts, easy to clean for healthier feeding

Our trusted Classic range has been improved to provide an even more enjoyable feeding experience. The clinically proven, anti-colic system is now integrated into the nipple, making the bottle extra easy to clean and assemble.

#### Clinically proven to reduce colic and discomfort\*

- Anti-colic system proven to reduce colic\*
- Fussing is reduced, especially at night

#### Easy to use, clean and assemble

- Leak-free for an enjoyable feeding
- Few parts for a quick and simple assembly
- Ergonomic shape for easy holding
- Easy to clean for healthier feeding

#### Other benefits

- BPA free
- Compatible with the Philips Avent range
- Different nipple flow rates available

# Highlights

## Less colic



Unlike other bottles, the clinically proven anti-colic system is now integrated into the nipple, making it easier to assemble the bottle correctly. As your baby feeds, the unique valve on the nipple opens to allow air into the bottle instead of your baby's tummy.

## Less fussing\*



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby bottle design affects "infant behavior". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min,  $p=0.05$ ). This was especially true during the night-time.\*

## Leak-free



The new Classic+ bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

## Easy to assemble in seconds



The new Classic+ bottle has only 4 parts for quick and simple assembly.

## Easy to hold



Due to the unique shape, the feeding bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

## BPA free


Following EU directive (2011/8/EU), the Philips Avent Classic+ feeding bottle is made of BPA free material (polypropylene).

## Compatibility across range



The Philips Avent Classic+ bottle is compatible with most of the Philips Avent range. We recommend to use the Classic+ bottle with Classic+ teats only.

## Fits your baby's growing needs



The Philips Avent Classic+ range offers different nipple flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All nipples are available in twin packs: Newborn, Slow, Medium, Fast and Variable Flow.

## Easy to clean



With only 4 parts, a wide bottle neck, and rounded corners, our bottle is easier to clean quickly and thoroughly. Ease of cleaning is important for a healthier feeding experience.

# Specifications

## Bottle

**Capacity:** 4oz/125ml

**Material:** BPA free\*

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Development stages

**Stages:** 0-12 months

## Ease of use

**Ease of use:** Dishwasher & microwave safe,  
Easy to assemble, Easy to clean, Easy to hold

## Functions

**Nipple:** Flexes to feeding rhythm, Easy latch on, Proven Anti-colic system

**Clinically proven:** clinically proven anti-colic

**Ease of use:** easy to clean and assemble,  
leak-free

## Material

**Bottle:** BPA free\*, Polypropylene

**Nipple:** BPA free\*, Silicone

## What is included

**Baby Bottle:** 4 pcs

