



Classic+ baby bottle

- 3 Bottles
- 9 oz/260 ml
- Slow flow teat
- 1m+

SCF573/35

Clinically proven to reduce colic and discomfort*

Designed for uninterrupted feeding

Our Classic+ bottle's Airflex venting system and textured teat are designed to minimise feeding interruptions and discomfort. With the integrated anti-colic valve, air is vented into the bottle and away from the baby's tummy.

Benefits

Clinically proven to reduce colic and discomfort*

- Anti-colic valve proven to reduce colic*
- 60% less fussing at night*
- Ribbed texture prevents collapse for uninterrupted feeding
- Teat shape designed for secure latch

Easy to use, clean and assemble

- Leak-free design
- Easy to clean and assemble with few parts

- Easy to hold
- Wide neck bottle with rounded corner for easy cleaning

Other benefits

- Compatible range from breastfeeding to cup
- This bottle is BPA-free
- Different teat flow rates available

PHILIPS


Features

Less colic



Our anti-colic valve is designed to keep air away from your baby's tummy to reduce colic and discomfort. As your baby feeds, the valve integrated into the teat flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce colic and discomfort.

Helps settle your baby



Philips Avent Anti-colic bottle reduces fussing. Babies fed with Philips Avent Anti-colic bottles experienced 60% less fussing at night than babies fed with a competitor's anti-colic bottle.*

Uninterrupted feeding



The teat shape allows a secure latch and the ribbed texture helps prevent teat collapse for uninterrupted, comfortable feeding.

Secure latch

Leak-free design

Our Anti-colic bottle is designed to prevent leakage whilst feeding for a truly enjoyable feeding experience.

Easy to assemble



Our Anti-colic bottle has few parts for quick and simple assembly.

Easy to hold



The unique bottle shape makes this bottle easy to hold and grip in any direction.

Easy to clean



With its wide neck and rounded corners, our bottle is easy to clean quickly and thoroughly. Ease of cleaning helps to ensure more hygienic feeding.

Growing up with your baby



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

BPA free



Philips Avent Anti-colic bottle is made of BPA free material (PP).

Fits your baby's growing needs



The Philips Avent Anti-colic bottle range offers different teat flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All teats are available in twin packs: Newborn, Slow, Medium, Fast and Variable Flow, and Thick feed.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Material

Bottle

BPA free
Polypropylene

Teat

BPA free
Silicone

What is included

Baby bottle 3 pcs

Bottle

Capacity 9 oz/260 ml

Design

Bottle design Ergonomic shape
Wide neck

Ease of use

Bottle use Dishwasher- and microwave-safe
Easy to assemble
Easy to clean
Easy to hold

Development stages

Stages 0-12 months

Functions

Anti-colic valve Airflex venting system designed to
reduce air ingestion
Teat Easy latch-on, Ribbed texture
prevents teat collapse, Proven anti-
colic system

Ease of use

easy to clean and assemble
leak-free design
4 pieces for easy assembly

Packaging dimensions

Height 21.00 cm
Width 16.00 cm
Depth 6.90 cm
Net weight 0.22 kg
Gross weight 0.26 kg
EAN 08710103699590
Number of products
included 1
Country of origin GB
Harmonised Systems
Code 392410

Outer Carton

Length 28.70 cm
Width 16.20 cm
Height 21.40 cm
Gross weight 1.15 kg
EAN 18710103699597
Number of consumer
packages 4

- * At 2 weeks of age, babies fed with a Philips Avent bottle showed less colic and significantly less fussing at night compared to babies fed with another competitor bottle.
- * Teat design proven to prevent teat collapse and associated air ingestion and feeding interruptions.
- * What is colic and how does it affect babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying and fussing.



data subject to change
2018, October 29
Version: 7.0.1
EAN: 08710103699590

© 2018 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com