

# PHILIPS

## AVENT

Avent Breast Milk  
Containers

Storage

SCF612/10



## Philips Avent storage system for easy storage Cup for milk storage

The Philips Avent storage system is a versatile, space saving storage system designed to grow with your baby. Use the same cup to store breast milk and baby food for your baby. It fits all Philips Avent breast pumps and nipples.

### Compatible with Philips Avent products

- Fit all Philips Avent breast pumps and nipples

### Easy to use, store and transport

- Easy to organize
- For use in fridge/freezer
- Leak-proof twist-on system
- Perfect for 'on-the-go'

# Specifications

## Country of origin

England

## What is included

Milk storage cup lid: 10 pcs

Milk storage cup adaptor: 2 pcs

Pre-sterilized VIA Cup (180 ml/ 6 oz): 10 pcs

## Material

BPA free\*

## Compact design

Ideal for going out

Ideal for travel

Stackable

## Development stages

Stages: 0 - 6 months, 6 - 12 months

## Functions

No leakage

Twist-on lid

