

PHILIPS

AVENT

Lids

SCF614/10



One system – many options

Storage system

VIA is a versatile, space saving storage system designed to grow with your baby. It allows you to use the same cup to store and feed baby food and breast milk to your baby. Fits all Philips AVENT breast pumps and nipples.

Benefits and features:

- Leak-proof, twist-on lid
- For use with the VIA Storage System

Specifications

Country of origin

England

What is included

Lid: 10 pcs

Material

BPA-Free

Compact design

Compact: Ideal for travel, Stackable

Development stages

Stages: 0-6 months, 0 - 6 months, 6 - 12 months

