

PHILIPS

AVENT

Avent Storage
System

Lids

SCF614/10



Philips Avent storage system for easy storage Cup for milk storage

The Philips Avent storage system is a versatile, space-saving storage system designed to grow with your baby. Use the same cup to store breast milk and baby food for your baby. It fits all Philips Avent breast pumps and teats.

Compatible with Philips Avent products

- Use with the Philips Avent breast milk containers

Easy to use, store and transport

- Leak-proof twist-on system
- Perfect for 'on-the-go'
- Store your baby's feed
- Take your baby's feed to go

Specifications

Country of origin

England

What is included

Milk storage cup lid: 10 pcs

Material

BPA free*

Compact design

Ideal for going out
Ideal for travel
Stackable

Development stages

Stages: 0–6 months, 0–6 months, 6–12 months

Functions

No leakage
Twist-on lid

