

PHILIPS

AVENT

Avent Refill Cups

180ml



SCF615/10



Philips Avent storage system for easy storage Cup for milk storage

The Philips Avent storage system is a versatile, space saving storage system designed to grow with your baby. Use the same cup to store breast milk and baby food for your baby. It fits all Philips Avent breast pumps and nipples.

Compatible with Philips Avent products

- Use with the Philip Avent breast milk containers

Easy to use, store and transport

- Leak-proof twist-on system
- Store your baby's feed
- Take baby's feed to go

Specifications

Country of origin

England

What is included

Pre-sterilized VIA Cup (180 ml/ 6 oz): 10 pcs

Material

BPA free*

Development stages

Stages: 0 - 6 months, 6 - 12 months

Convenience

Travel ease: Compact

Functions

No leakage

Easy storage: Store/freeze milk and food

Twist-on lid

