

# PHILIPS

## AVENT

### Refill Cups

240ml



SCF616/05



## One system – many options

### Storage system

VIA is a versatile, space saving storage system designed to grow with your baby. It allows you to use the same cup to store and feed baby food and breast milk to your baby. Fits all Philips AVENT breast pumps and nipples.

#### **Benefits and features:**

- Easy to organize
- Leak-proof, twist-on lid
- For use with the VIA Storage System
- Perfect for 'on-the-go'
- For use in fridge/freezer
- Fits all Philips AVENT breast pumps and nipples
- Pre-sterilized VIA Cups are ready to use

# Specifications

**Country of origin**

England

**What is included**

Storage Cup (240 ml/8oz): 5 pcs

**Development stages**

Stages: 0-6 months, 0 - 6 months, 6 - 12 months

**Material**

BPA-Free

