

# PHILIPS

## AVENT

### Breast milk storage cups

Includes convenient adaptors

180ml/6oz

10 pcs



SCF618/10



## Store breast milk securely

### With leak proof lid

Express, store and feed your breast milk efficiently, with our new storage cup. Sterilize and re-use the storage cup with the Philips Avent pump or teats-one system, many options!

#### **Leak proof lid**

- For safe storage and transport

#### **Write-on cups and lids**

- To track easy dates & contents

#### **Stackable cups and lids**

- To track easily dates and contents
- An organized fridge & freezer

#### **Maximal convenience**

- Convenient adaptors included
- One system- many options
- For use in fridge and freezer
- Easy to use and clean
- This product is BPA Free

# Highlights

## Leak proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

## Write-on cups and lids



To track easily dates and contents.

## For easy storage

For easy storage.

## Stackable cups and lids



For an organized fridge and freezer.

## Convenient adaptors included



For expressing and feeding.

## One system- many options



Compatible with Philips Avent breast pumps, Classic, Classic + and Natural teats.

## For use in fridge and freezer



For maximum flexibility.

## Easy to use and clean



Can be used safely in bottle warmer, microwave, dishwasher and sterilizer.

## BPA Free

This product is BPA Free

# Specifications

## Weight and dimensions

Retail pack dimensions: 158x194x80 mm

## Country of origin

Designed in: UK

Produced in: UK

## What is included

180ml/6oz storage cups: 10 (re-usable) pcs

Lids: 10 pcs

Adapter: 2x

## Development stages

Stages: 0+ months

