

PHILIPS

AVENT

Breast milk storage cups

Includes convenient adaptors

180ml/6oz

10 pcs



SCF618/10



Store breast milk securely

With leak proof lid

Express, store and feed your breast milk efficiently, with our new storage cup. Sterilize and re-use the storage cup with the Philips Avent pump or teats-one system, many options!

Leak proof lid

- For safe storage and transport

Write-on cups and lids

- To track easy dates & contents

Stackable cups and lids

- To track easily dates and contents
- An organized fridge & freezer

Maximal convenience

- Convenient adaptors included
- One system- many options
- For use in fridge and freezer
- Easy to use and clean
- This product is BPA Free

Highlights

Leak proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Write-on cups and lids



To track easily dates and contents.

For easy storage

For easy storage.

Stackable cups and lids



For an organized fridge and freezer.

Convenient adaptors included



For expressing and feeding.

One system- many options



Compatible with Philips Avent breast pumps, Classic, Classic + and Natural teats.

For use in fridge and freezer



For maximum flexibility.

Easy to use and clean



Can be used safely in bottle warmer, microwave, dishwasher and sterilizer.

BPA Free

This product is BPA Free

Specifications

Weight and dimensions

Retail pack dimensions: 158x194x80 mm

Country of origin

Designed in: UK

Produced in: UK

What is included

180ml/6oz storage cups: 10 (re-usable) pcs

Lids: 10 pcs

Adapter: 2x

Development stages

Stages: 0+ months

