



Breast milk storage cup

- 180ml/6oz
- 5 pcs

SCF619/05

Store breast milk securely

With leak proof lid

Express, store and feed your breast milk efficiently, with our new storage cup. Sterilize and re-use the storage cup with the Philips Avent pump or teats- one system, many options!

Benefits

Leak proof lid

- For safe storage and transport

Write-on cups and lids

- To track easy dates & contents

Stackable cups and lids

- To track easily dates and contents
- An organized fridge & freezer

Maximal convenience

- One system- many options
- For use in fridge and freezer
- Easy to use and clean
- This product is BPA Free

PHILIPS

Features

Leak proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Write-on cups and lids



To track easily dates and contents.

For easy storage

For easy storage.

Stackable cups and lids



For an organized fridge and freezer.

One system- many options



Compatible with Philips Avent breast pumps, Classic, Classic + and Natural teats.

For use in fridge and freezer



For maximum flexibility.

Easy to use and clean



Can be used safely in bottle warmer, microwave, dishwasher and sterilizer.

BPA Free

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Country of origin

Designed in	UK
Produced in	UK

Development stages

Stages	0+ months
--------	-----------

What is included

180ml/6oz storage cups	5 (re-usable) pcs
Lids	5 pcs



data subject to change
2019, January 10
Version: 4.2.1
UPC: 00075020041289

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com