

# PHILIPS

# AVENT

Natural baby bottle

1 Bottle

9oz/260ml

Slow flow teat

Natural teat shape



SCF627/22



## The most natural way to bottle feed

### Easy to combine with breastfeeding

Our new bottle helps to make bottle feeding more natural for your baby and you. The teat features an innovative petal design for natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

#### **Easy to combine breast and bottle feeding**

- Natural latch on due to the wide breast shaped teat
- This bottle is BPA free\*
- Trusted by moms since 1984

#### **A more comfortable and contented feed for baby**

- Unique petals for a soft, flexible teat without collapse

#### **Vents air away from your baby's tummy**

- Advanced anti-colic system with innovative twin valve

#### **Easy to hold and grip in any direction**

- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly

# Highlights

## Natural latch on



The wide breast shaped teat promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.

## Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

## This bottle is BPA free\*



The Philips Avent Natural bottle is made of BPA free\* material (polypropylene).

## Ergonomic shape



Extra comfortable shape, easy to hold and grip in any direction, even for baby's tiny hands.

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

## Trusted by moms

Since 1984, Philips Avent has been designing and manufacturing products that are inspired by nature and have been developed through extensive research and clinical trials.

# Specifications

## Country of origin

Indonesia

## Material

Bottle: Polypropylene, BPA free\*

Teat: BPA free\*, Silicone

## What is included

Baby Bottle: 1 pcs

## Bottle

Material: BPA free\*

Capacity: 9oz/260ml

## Design

Bottle design: Ergonomic shape, Wide neck

## Ease of use

Ease of use: Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

## Development stages

Stages: 0-12 months

## Functions

Latch on: Easy combine breast and bottle

Teat: Advanced anti-colic system, Extra soft and flexible teat, Natural latch on, Unique comfort petals

