

PHILIPS

AVENT

Classic+ nipple

2 pieces

Variable flow

3m+

SCF635/02



Clinically proven to reduce colic and discomfort*

Trusted for 30 years

With the Philips Avent Classic+ nipple, your baby controls milk flow, which can help reduce over eating and spit-up. The unique valve flexes to allow air into the bottle, promoting healthy, active feeding and reducing fussing and colic.*

Ideal for thicker liquids - 3M+

- Adjust the flow rate to baby's convenience

Vents air away from your baby's tummy

- Unique anti-colic system

0% BPA (Bisphenol A Free)

- BPA free nipple

More comfortable feed at baby's chosen pace

- Easy latch on due to the unique valve on the nipple

Reduce fussing and discomfort

- Clinically proven to significantly reduce fussing

Highlights

Easy latch on



Unique valve on the nipple flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby to minimize overeating and spit-up, burping and gas

Significantly reduces fussing



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby bottle design affects "infant behavior". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min, $p=0.05$) This was especially true during the night-time.**

Unique anti-colic valve



As your baby feeds, the unique valve on the nipple flexes to allow air into the bottle instead of your baby's tummy.

BPA free



This nipple is made from silicone - a BPA-free material (Following EU directive 2011/8/EU)

Fit your baby's growing needs




The Philips Avent Classic bottle offers different flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All nipples are available in twin packs.

Compatible with Classic bottle

We advise to use the Classic nipple with Classic baby bottles only.

Adjust the flow rate



With the variable flow nipple, you can adjust the flow rate according to liquid's thickness and suit baby's feeding rhythm perfectly. We recommend to use this nipple for thicker liquids such as AR (thickened) milk, milk mixed with baby rice, milk mixed with baby food, soup... Once baby is ready for even thicker food, we recommend to use the Philips Avent Thick Feed nipple.

Specifications

Country of origin

England

Material

Nipple: Silicone, BPA free*

What is included

Soft, Variable Flow Nipple: 2 pcs

Design

Nipple design: Wide

Nipple

Flow speed: Variable flow

Holes: Slot hole

Months: 3m+ (vari-flow)

Development stages

Stage: 0-6 months

Functions

Nipple: Flexes to feeding rhythm

Latch on: Easy latch on

Anti-colic valve: Two piece anti-colic system



* 0% BPA, following EU regulation 10/2011

** A clinical study showed that at two weeks of age, babies showed less fussing than babies fed with another bottle. (Study conducted by the Institute of Child Health, London, 2008.)

*** A clinical study showed that at 2 weeks of age, infants fed with an Avent bottle showed less colic than infants fed with another leading bottle, especially at night.