

PHILIPS

AVENT

Food storage cup

Standard



SCF639/05



Ideal food storage for home and away With leak proof lid

Store and feed your healthy weaning meals with our new food storage cup. Sterilize and re-use the storage cup with the Philips Avent at home and away.

Leak proof lid

- For safe storage and transport

Write-on cups and lids

- To track easy dates & contents

Stackable cups and lids

- To track easily dates and contents
- An organized fridge & freezer

Maximal convenience

- For use in fridge and freezer
- Easy to use and clean
- This product is BPA Free

Highlights

Leak proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Write-on cups and lids



To track easily dates and contents.

For easy storage



For easy storage.

Stackable cups and lids



For an organized fridge and freezer.

For use in fridge and freezer



For maximum flexibility.

Easy to use and clean



Can be used safely in bottle warmer, microwave, dishwasher and sterilizer.

BPA Free

This product is BPA Free

Specifications

Country of origin

Designed in: UK
Produced in: UK

What is included

Lids: 5 pcs
240ml/8oz storage cups: 5 (re-usable) pcs

Development stages

Stages: 6 months +

Weight and dimensions

Retail pack dimensions: 80x80x200mm
(WxDxH)

