

PHILIPS

AVENT

Natural nipple

2 pieces

Newborn flow

0m+



SCF651/27



The natural way to bottle feed

Avent nipple with petal design

Our new nipple helps to make bottle feeding more natural for your baby and you. The nipple features an innovative petal design for natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

0% BPA (Bisphenol A Free)

- BPA free nipple

Other benefits

- Different flow rates for the most comfortable feed
- Compatible with Philips Avent Natural feeding bottle



Highlights

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

BPA free



This nipple is made from silicone - a BPA-free material (Following EU directive 2011/8/EU)

Different flow rates

Philips Avent offers four different flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. Philips Avent offers the newborn flow nipple with 1 hole (0M+), the slow flow nipple with 2 holes (1M+), the medium flow nipple with 3 holes (3M+) and the fast flow nipple with 4 holes (6M+). All nipples are available in twin packs.

Compatible with Natural bottle

We advise to use the Natural feeding nipples with Natural bottles only.

Specifications

Country of origin

England

Material

Nipple: Silicone, BPA free*

What is included

Soft, Newborn Flow nipple: 2 pcs

Design

Nipple design: Breast shaped nipple, Unique comfort petals, Extra wide

Nipple

Flow speed: Newborn flow

Holes: 1 hole

Months: 0m+

Development stages

Stages: 0-6 months

Functions

Latch on: Easy combine breast and bottle, Natural latch on

Nipple: Unique comfort petals, Extra soft and flexible teat

Anti-colic valve: Advanced anti-colic system

