

PHILIPS

AVENT

Natural nipple

2 pieces

Medium flow

3m+



SCF653



The most natural way to bottle feed

Natural latch on

Our Natural ultra soft nipple more closely resembles the breast. The wide breast-shaped nipple with flexible spiral design and comfort petals allows natural latch on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft nipple designed to mimic the feel of the breast

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Unique anti-colic Airflex vent technology

Other benefits

- Nipples with different flow rates available

Highlights

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Anti-colic vent technology



Our unique anti-colic Airflex vent technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Choosing the right nipple

The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Specifications

Material

Nipple: Silicone, BPA free*

What is included

Ultra soft Medium flow nipple: 2 pcs

Development stages

Stages: 3m+

Functions

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Ultra soft and flexible nipple

Anti-colic valve: Advanced anti-colic system

