

PHILIPS

AVENT

Natural nipple

2 pieces

Variable flow

3m+

SCF655



The most natural way to bottle feed

Natural latch on

With the Natural variable flow nipple, you can adjust the flow rate by simply turning the bottle. The comfort petals and natural nipple shape allow natural latch on and make it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft and smooth silicone for your baby's changing needs

A comfortable and contented feed for your baby

- Flexible, anti-collapse ribbed teat design

Vents air away from your baby's tummy

- Unique anti-colic Airflex vent technology

Other benefits

- Nipple with different flow rates available

Ideal for thicker liquids - 3M+

- Adjustable flow rate to suit baby's feeding rhythm

Highlights

Specifications

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Soft and smooth silicone

Bite resistant and smooth silicone nipple designed for your baby's changing needs.

Flexible, anti-collapse design



Petals and ridges inside the nipple increase softness and flexibility without nipple collapse.

Anti-colic vent technology



Our unique anti-colic Airflex vent technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Choosing the right nipple

The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Adjust the flow rate

With the variable flow nipple, you can adjust the flow rate according to the liquid's thickness and suit the baby's feeding rhythm perfectly.

We recommend to use this nipple for thicker liquids such as AR (thickened) milk, milk mixed with baby rice, milk mixed with baby food, soup...

Material

Nipple: Silicone, BPA free*

What is included

Soft Variable Flow nipple: 2 pcs

Development stages

Stage: 3m+

Functions

Anti-colic valve: Advanced anti-colic system

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Soft and anti-collapse

