

PHILIPS

AVENT

Feeding Bottles

Advanced Classic

11 oz



SCF666/17



Good feeding, good sleeping

Helps settle your baby, especially at night**

Fussing is by far the most commonly known crying behaviour in babies. The Philips AVENT Advanced Classic feeding bottle significantly reduces colic*** and fussing. Fussing is reduced especially at night.**

Reduce fussing and discomfort

- Clinically proven to significantly reduce fussing

One-piece anti-colic valve system

- Allows air into the bottle instead of your baby's tummy

Highly durable

- Honey-coloured BPA-Free PES material for extra durability

Easy to clean

- Only three components for ease of cleaning

To fit your baby's growing needs

- Five different teat flow rates are available

Highlights

Helps settle your baby




Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant feeding bottle design affects "infant behaviour". The Philips AVENT Classic baby feeding bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs. 74 min, $p=0.05$). This was especially true at night.**

Unique anti-colic valve



As your baby feeds, the unique skirt on the AVENT Teat flexes to allow air into the bottle instead of your baby's tummy. Like breastfeeding, your baby controls the milk flow.

Fits your baby's growing needs



Five different teat flow rates are available.

BPA Free (0% Bisphenol A)



This bottle is made from PES - a BPA-Free material

Easy to clean

The Philips AVENT feeding bottle consists of only three components for ease of cleaning and has a dome cap for hygienic storage and transportation. In a randomised clinical trial, the design of the Philips AVENT bottle was compared with another leading brand over a 12-month period. The maternal satisfaction in terms of ease of cleaning ($p=0.04$) and ease of assembly ($p=0.02$) was significantly better for the Philips AVENT bottle.

Honey-coloured PES material

This BPA-Free bottle is made of honey-coloured PES material for extra durability

Specifications

Country of origin

England

Material

BPA free*

What is included

Advanced Classic Bottle: 1 pcs
Classic variable flow teat: 1 pcs

Bottle

Capacity: 330 oz

Material: Highly durable, BPA free

Design

Bottle design: Easy to clean, Easy to assemble, Easy to hold, Wide neck

Functions

Anti-colic valve: One piece valve

Development stages

Stages: 0–6 months



* A clinical study showed that at two weeks of age, babies showed less fussing than babies fed with another bottle. (Study conducted by the Institute of Child Health, London, 2008.)

** A clinical study showed that at 2 weeks of age, infants fed with an AVENT bottle showed less colic than infants fed with another leading bottle, especially at night.