

PHILIPS

AVENT

Feeding bottle

Advanced Classic

11oz

Medium Flow Nipple



SCF666/27



Good feeding, good sleeping

Helps settle your baby, especially at night**

Fussing is by far the most commonly known crying behavior in babies. The Philips AVENT Advanced Classic feeding bottle significantly reduces colic*** and fussing. Fussing is reduced especially at night.**

Reduce fussing and discomfort

- Clinically proven to significantly reduce fussing**

One piece anti-colic valve system

- Allows air into the bottle instead of your baby's tummy

Highly durable

- Honey-colored BPA-Free PES material for extra durability

Easy to clean

- Only three components for ease of cleaning

To fit your baby's growing needs

- Five different nipple flow rates are available

Highlights

Significantly reduces fussing



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant feeding bottle design affects "infant behaviour". The Philips AVENT Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min, $p=0.05$). This was especially true during night-time.

Unique anti-colic valve




As your baby feeds, the unique skirt on the AVENT Nipple flexes to allow air into the bottle instead of your baby's tummy. Like natural breastfeeding your baby controls the milk flow.

0% BPA (Bisphenol A Free)



This bottle is made from PES - a BPA-Free material

Fits your baby's growing needs



Five different nipple flow rates are available.

Easy to clean

The Philips AVENT feeding bottle consists of only three components for ease of cleaning and has an dormal cap for hygienic storage and transportation. In a randomised clinical trial the design of the Philips AVENT bottle was compared with another leading brand over a 12 month period. The maternal satisfaction in terms of ease of cleaning ($p=0.004$) and ease of assembly ($p=0.02$) was significantly better in the Philips AVENT bottle.

Honey-colored PES material

This BPA-Free bottle (0% Bisphenol A) is made of honey-colored PES material for extra durability

Specifications

Country of origin

England

Material

BPA free*

What is included

Advanced Classic Bottle: 2 pcs

Bottle

Capacity: 330 ml

Material: BPA free, Highly durable

Design

Bottle design: Easy to assemble, Easy to clean, Easy to hold, Wide neck

Development stages

Stages: 0 - 6 months

Functions

Anti-colic valve: One piece valve



** A clinical study showed that at two weeks of age, babies showed less fussing than babies fed with another bottle. (Study conducted by the Institute of Child Health, London, 2008.)

*** A clinical study showed that at 2 weeks of age, infants fed with an AVENT bottle showed less fussing than infants fed with another leading bottle, especially at night.