

PHILIPS

AVENT

Natural glass baby bottle

3 Bottles

4oz/120ml

Newborn flow nipple

0m+



SCF671



Easy to combine with breastfeeding

Natural latch on

Our Natural glass bottle helps to make bottle feeding more natural for baby and you. The nipple features a comfort petal design for natural latch on similar to the breast, making it easy for baby to combine breast and bottle feeding.

Pure and safe

- Heat resistance
- Pharmaceutical grade glass

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse
- Compatible with Philips Avent range

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

Other benefits


- Simple to use and clean, quick and easy assembly
- This bottle is BPA free*



asimpleswitch.com

Highlights

Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

Pharmaceutical grade glass


Premium quality borosilicate glass to ensure the highest quality for ultimate purity.

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

Growing with your baby

The Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (glass, polypropylene and silicone).

Specifications

Material

Bottle: BPA free*, Glass
Nipple: BPA free*, Silicone

What is included

Baby Bottle: 3 pcs

Bottle

Material: BPA free*
Capacity: 4oz/120ml

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Ease of use: Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

Development stages

Stages: 0-12 months

Functions

Latch on: Easy combine breast and bottle
Nipple: Unique comfort petals, Advanced anti-colic system, Extra soft and flexible nipple, Natural latch on



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

