

# PHILIPS

## AVENT

Natural glass baby bottle

1 Bottle

8oz/240ml

Slow flow nipple

1m+



SCF673/17



## Easy to combine with breastfeeding Natural latch on

Our Natural glass bottle helps to make bottle feeding more natural for baby and you. The nipple features a comfort petal design for natural latch on similar to the breast, making it easy for baby to combine breast and bottle feeding.

### Pure and safe

- Heat resistance
- Pharmaceutical grade glass

### Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple

### A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse
- Compatible with Philips Avent range

### Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

### Other benefits

- Simple to use and clean, quick and easy assembly
- This bottle is BPA free\*



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

## Pharmaceutical grade glass


Fully recyclable, premium quality borosilicate glass to ensure the highest quality for ultimate purity.

## Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

## Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

## Growing with your baby

The Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

## This bottle is BPA free\*



The Philips Avent Natural bottle is made of BPA free\* material (glass, polypropylene and silicone).

# Specifications

## Country of origin

Assembled in UK or in the USA: With components from the USA and UK

## Material

**Bottle:** BPA free\*, Glass  
**Nipple:** Silicone, BPA free\*  
**BPA free\*:** Bottle made of Borosilicate type 1 glass; Nipple made of silicone; Screw ring made of PP (Polypropylene)

## What is included

**Baby Bottle:** 1 pcs

## Bottle

**Material:** BPA free\*  
**Capacity:** 8oz/240ml

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Ease of use:** Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

## Development stages

**Stages:** 0-12 months

## Functions

**Latch on:** Easy combine breast and bottle  
**Nipple:** Unique comfort petals, Advanced anti-colic system, Extra soft and flexible nipple, Natural latch on

## Dimensions & Weight

**Dimensions:** Height: 8.7cm ; Diametre: 6.6cm  
**Weight:** 1 bottle: 14.67g



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

\* 0% BPA, following EU regulation 10/2011