

PHILIPS

AVENT

Natural glass baby bottle

2 Natural Glass bottles

8oz

Slow Flow Nipple



SCF673/27



The natural way to bottle feed

Avent nipple with petal design

Our Natural Glass bottle helps to make bottle feeding more natural for baby and you. The nipple features an innovative petal design for natural latch on similar to the breast, making it easy for baby to combine breast and bottle feeding.

Pure and safe

- Heat resistance
- Pharmaceutical grade glass

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

Other benefits

- Simple to use and clean, quick and easy assembly
- This bottle is BPA free*



asimpleswitch.com

Highlights

Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

Pharmaceutical grade glass

Premium quality borosilicate glass to ensure the highest quality for ultimate purity.

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (glass, polypropylene and silicone).

Specifications

Bottle

Material: BPA free*
Capacity: 8oz/240ml

What is included

Baby Bottle: 2 pcs

Material

BPA free*: Bottle made of Borosilicate type 1 glass; Nipple made of silicone; Screw ring made of PP (Polypropylene)

Bottle: BPA free*, Glass

Nipple: BPA free*, Silicone

Dimensions & Weight

Dimensions: Height: 8.7cm ; Diametre: 6.6cm

Weight: 1 bottle: 14.67g

Functions

Latch on: Easy combine breast and bottle

Nipple: Unique comfort petals, Advanced anti-colic system, Extra soft and flexible nipple, Natural latch on

Ease of use

Ease of use: Dishwasher & microwave safe, Easy to assemble, Easy to clean, Easy to hold

Design

Bottle design: Ergonomic shape, Wide neck

Development stages

Stage: 0-6 months



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.