

PHILIPS

AVENT

Avent Breast Milk
Containers

4 oz x 4



SCF680/04



Easy breast milk storage and feeding No need to transfer milk

Philips Avent breast milk container SCF680/04 for storage in the fridge or freezer. Express directly into the containers by attaching them to your breast pump. Simply replace the sealing discs with a teat for feeding. (Teats not included)

Express

- Express milk directly into storage container

Store

- Insert a sealing disc to store breast milk

Feed

- Simply replace the sealing disc with an Avent teat

BPA Free (0% Bisphenol A)

- This bottle is made from PP - a BPA-Free material

The Philips Avent Advantage

- All Avent teat and spouts can be used

Highlights

The Philips Avent Advantage



All Philips Avent teats and spouts can be used with the Philips Avent Breast Milk Container. (Teats not included.)

BPA Free (0% Bisphenol A)



This bottle is made from PP - a BPA-Free material

Express



Express milk directly into storage container using any Philips Avent Breast Pump

Store



Simply insert a sealing disc to store breast milk in the fridge or freezer

Feed



Replace the sealing disc with an Avent teat for feeding - no need to transfer milk. (Teat is not included.)

Specifications

What is included

Breast milk container (125 ml/4 oz): 4 pcs
Sealing disc for milk storage: 4 pcs
Screw ring: 4 pcs

Design

Bottle design: Easy to assemble, Easy to clean, Easy to hold, Wide neck
Ease of use: Compact

Storage system

Storage: Breast milk storage, Dishwasher safe, Freezer safe, Sterilisable, Use in fridge, Works with Philips Avent range

Development stages

Stages: 0-6 months

Material

BPA free*

