

PHILIPS

AVENT

Baby Bottle Trainer
Kit

9oz

Variable Flow Nipple

SCF683/07



Baby's first step to using a cup

Easy-grip handles for bottles

Designed to be gentle on baby's gums, the Philips AVENT bottle trainer SCF683/07 with a soft, non-spill drinking spout

Help your baby learn to drink independently

- Attach the easy-grip handles to the feeding bottle

Gentle on baby's gums

- Soft non-spill drinking spout

Easy to grip handles

- Easy to grip, removable handles

0% BPA (Bisphenol A Free)

- This bottle is made from PP - a BPA Free material

Transform into a trainer cup

- Replace nipple with the soft spout

Highlights

Attach the easy grip handles



Attach the easy-grip handles to the Philips AVENT baby bottle.

Easy to grip, removable handles



Attach the easy-grip handles to the baby bottle to help your baby learn to drink independently

Soft non-spill drinking spout



Replace the nipple with the soft, non-spill drinking spout - designed to be gentle on your baby's gums

Replace with soft spout



Transform feeding bottle by simply replacing the nipple with the soft, non-spill drinking spout

0% BPA (Bisphenol A Free)



This bottle is made from PP - a BPA Free material

Specifications

Development stages

Stages: 0 - 6 months, 6 months +

Material

BPA free*

Bottle

Capacity: 260 oz

Material: BPA free

Design

Bottle design: Easy to assemble, Easy to clean, Easy to hold, Wide neck

What is included

Natural Feeding Bottle: 1 pcs

Easy grip handles: 1 pcs

Soft, Variable Flow Nipple: 1 pcs

Soft spout: 1 pcs

Dome cap: 1 pcs

Screw ring: 1 pcs

