

PHILIPS

AVENIT

Natural baby bottle

1 Bottle

4oz/125ml

Newborn flow nipple

0m+



SCF690/17



Easy to combine with breastfeeding Natural latch on

Our new bottle helps to make bottle feeding more natural for your baby and you. The nipple features an innovative petal design for natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.



Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Trusted by moms since 1984
- This bottle is BPA free*

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse
- Compatible with Philips Avent range
- Available in different sizes (4oz/125ml and 9oz/260ml)

Air goes into the bottle, not baby's tummy

- Advanced anti colic system with innovative twin valve

Easy to hold and grip in any direction

- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly



Highlights

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve design reduces colic and discomfort by venting air into the bottle and not baby's tummy.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Trusted by moms

Since 1984, Philips Avent has been designing and manufacturing products that are inspired by nature and have been developed through extensive research and clinical trials.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Growing with your baby

The Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

Different sizes

The Philips Avent Natural bottle is available in 3 sizes: 4oz/125ml, 9oz/260ml and 11oz/330ml. All bottles are available in single and multi packs

This bottle is BPA free*



The new Philips Avent Natural bottle is made of BPA free* material (polypropylene).



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO₂ emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Bottle

Capacity: 4oz/125ml

Material: BPA free*

Country of origin

England

Design

Bottle design: Ergonomic shape, Wide neck

Development stages

Stages: 0-12 months

Ease of use

Ease of use: Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

Functions

Nipple: Advanced anti-colic system, Natural latch on, Extra soft and flexible nipple, Unique comfort petals

Latch on: Easy combine breast and bottle

Ease of use: easy to clean and assemble

Material

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

What is included

Baby Bottle: 1 pcs

