

# PHILIPS

## AVENT

### Natural baby bottle

2 Bottles

4oz/125ml

Newborn flow nipple

0m+



SCF690/27



## Easy to combine with breastfeeding Natural latch on

Our new bottle helps to make bottle feeding more natural for your baby and you. The nipple features an innovative petal design for natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

#### Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Trusted by moms since 1984
- This bottle is BPA free\*

#### A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse
- Compatible with Philips Avent range
- Available in different sizes

#### Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

#### Easy to hold and grip in any direction

- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

## Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

## Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

## Trusted by moms

Since 1984, Philips Avent has been designing and manufacturing products that are inspired by nature and have been developed through extensive research and clinical trials.

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

## Growing with your baby

The Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

## Different sizes

The Philips Avent Natural bottle is available in 4 sizes: 2oz/60ml, 4oz/125ml, 9oz/260ml and 11oz/330ml. Bottles are available in single and multi packs.

## This bottle is BPA free\*



The Philips Avent Natural bottle is made of BPA free\* material (polypropylene).



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Bottle

**Capacity:** 4oz/125ml

**Material:** BPA free\*

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Development stages

**Stages:** 0-12 months

## Ease of use

**Ease of use:** Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

## Functions

**Nipple:** Unique comfort petals, Advanced anti-colic system, Extra soft and flexible nipple, Natural latch on

**Latch on:** Easy combine breast and bottle

## Material

**Bottle:** Polypropylene, BPA free\*

**Nipple:** Silicone, BPA free\*

## What is included

**Baby Bottle:** 2 pcs

