

PHILIPS

AVENT

Natural baby bottle

2 Bottles

11oz/330ml

Medium flow nipple

3m+



SCF696/27



Easy to combine with breastfeeding

Natural latch on

Our new bottle helps to make bottle feeding more natural for your baby and you. The nipple features an innovative petal design for natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.



Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- This bottle is BPA free*
- Trusted by moms since 1984

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse
- Available in different sizes
- Compatible with Philips Avent range

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

Easy to hold and grip in any direction

- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly



asimpleswitch.com

Highlights

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).

Different sizes

The Philips Avent Natural bottle is available in 4 sizes: 2oz/60ml, 4oz/125ml, 9oz/260ml and 11oz/330ml. Bottles are available in single and multi packs.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Growing with your baby

The Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

Trusted by moms

Since 1984, Philips Avent has been designing and manufacturing products that are inspired by nature and have been developed through extensive research and clinical trials.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Material

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Baby Bottle: 2 pcs

Bottle

Material: BPA free*

Capacity: 11oz/330ml

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Ease of use: Easy to assemble, Easy to clean,
Easy to hold, Dishwasher & microwave safe

Development stages

Stage: 0-12 months

Functions

Latch on: Easy combine breast and bottle

Nipple: Unique comfort petals, Advanced anti-colic system, Extra soft and flexible nipple, Natural latch on

