

PHILIPS

AVENT

Baby Bottle

1 Natural bottle

2 oz/60 ml

First flow teat

Natural teat shape



SCF699/17



Ideal for newborns

Easy to combine with breastfeeding

Our Natural 2-oz bottle provides the ideal amount of food for newborns. Designed with the First Flow teat, it has a slower, more controlled flow rate. The soft, textured, wide-shaped teat more closely resembles the breast.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast-shaped teat
- This bottle is BPA free*

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible teat without collapse
- Available in different sizes
- Compatible with Philips Avent range

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

Ideal for newborns

- Slower-flow teat
- smaller 2-oz/60-ml bottle designed for newborn feeding
- Soft, textured teat



asimpleswitch.com

Highlights

Slower-flow teat



A smaller hole size delivers a more controlled flow rate for babies who drink more slowly. The Natural First teat provides the ideal flow rate for starting to bottle feed your baby.

2-oz/60-ml bottle container



A smaller-size bottle to help ensure the appropriate amount for baby's smaller tummy

Soft, textured teat



The softer, textured teat more closely resembles the breast.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

Different sizes

The Philips Avent Natural bottle is available in 4 sizes: 2 oz/60 ml, 4 oz/125 ml, 9 oz/260 ml and 11 oz/330 ml. Bottles are available in single and multi packs.

Growing with your baby

The Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise using the Natural bottles with Natural feeding teats only.

Natural latch on




The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA-free* material (polypropylene).

Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Material

Bottle: Polypropylene, BPA-free*

Teat: BPA-free*, Silicone

BPA-free*

What is included

Baby Bottle: 1 pcs

Bottle

Material: BPA-free*

Capacity: 2 oz/60 ml

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Ease of use: Easy to assemble, Easy to clean,
Dishwasher- and microwave-safe

Development stages

Stage: 0-6 months

Functions

Latch on: Easily combine breast and bottle

Teat: Unique comfort petals, Advanced anti-colic system, Extra soft and flexible teat,
Natural latch on

