

PHILIPS

AVENT

Natural glass baby bottle

3 Bottles

4oz/120ml

Newborn flow nipple

0m+



SCF701/37



The most natural way to bottle feed

Natural latch on

Our Natural bottle with an ultra soft nipple more closely resembles the breast. The wide breast-shaped nipple with flexible spiral design and comfort petals allows natural latch on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Soft nipple designed to mimic the feel of the breast
- Natural latch on due to the wide breast shaped nipple
- Flexible spiral design combined with comfort petals

Pure glass

- Heat resistance
- Pharmaceutical grade glass

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- Simple to use and clean, quick and easy assembly
- Nipples with different flow rates available
- Compatible with the Philips Avent Natural range
- This bottle is BPA free*



asimpleswitch.com

Highlights

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

Pharmaceutical grade glass

Premium quality borosilicate glass to ensure the highest quality for ultimate purity.

Anti-colic valve technology

Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Choosing the right nipple



The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Development stages

Stages: 0-6 months

Ease of use

Bottle use: Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

Functions

Anti-colic valve: Advanced anti-colic system

Latch on: Easy combine breast and bottle, Natural latch on

Nipple: Unique comfort petals, Ultra soft and flexible nipple

Material

Bottle: BPA free*, Glass

Nipple: BPA free*, Silicone

What is included

Baby bottle: 3 pcs

