

PHILIPS

AVENT

Food storage cup

Set



SCF721/20



Ideal food storage for home and away With leak proof lid

Philips Avent food storage cups SCF720/10 are versatile, easy to use and designed to grow with your baby. The cups are ideal for storing food both at home and away.

Leak proof lid

- For safe storage and transport

Write-on cups and lids

- To track easy dates & contents

Stackable cups and lids

- To track easily dates and contents
- An organized fridge & freezer

Maximal convenience

- For use in fridge and freezer
- Easy to use and clean
- Includes weaning spoon
- This product is BPA Free

Highlights

Specifications

Leak proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

For use in fridge and freezer



For maximum flexibility.

Country of origin

Designed in: UK
Produced in: UK

What is included

180ml/6oz storage cups: 10 (re-usable) pcs
240ml/9oz storage cups: 10 (re-usable) pcs
Lids: 20 pcs
Weaning spoon: 1 pcs

Development stages

Stages: 6 months +

Write-on cups and lids



To track easily dates and contents.

Easy to use and clean



Can be used safely in bottle warmer, microwave, dishwasher and sterilizer.

For easy storage



For easy storage.

Includes weaning spoon



To wean your baby home made meals.

Stackable cups and lids



For an organized fridge and freezer.

BPA Free

This product is BPA Free

