

# PHILIPS

## AVENT

Anti-colic baby  
bottle

2 Bottles

11oz/330ml

Medium flow teat

3m+

SCF816/27



## Clinically proven to reduce colic and discomfort\*

### Designed for uninterrupted feeding

Our Anti-colic bottle is designed to keep air away from your baby's tummy to reduce colic and discomfort, while also minimizing feeding interruptions. With its integrated anti-colic valve, air is vented away from the baby's tummy.

#### Clinically proven to reduce colic and discomfort\*

- Anti-colic valve proven to reduce colic\*
- 60% less fussing at night\*
- Ribbed texture prevents collapse for uninterrupted feeding
- Teat shape designed for secure latch

#### Easy to use, clean and assemble

- Easy to clean and assemble with few parts
- Easy to hold
- Wide neck bottle with rounded corner for easy cleaning
- Leak-free design

#### Other benefits

- Compatible range from breastfeeding to cup
- This bottle is BPA free
- Different teat flow rates available

# Highlights

## Less colic



Our anti-colic valve is designed to keep air away from your baby's tummy to reduce colic and discomfort. As your baby feeds, the valve integrated into the teat flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce colic and discomfort.

## Less fussing



Philips Avent Anti-colic bottle reduces fussing. Babies fed with Philips Avent Anti-colic bottles experienced 60% less fussing at night, than babies fed with a competitor's anti-colic bottle.\*

## Uninterrupted feeding



The teat shape allows secure latch and the ribbed texture helps prevent teat collapse for uninterrupted, comfortable feeding.

## Secure latch

Teat shape designed for secure latch

## Easy to assemble



Our Anti-colic bottle has few parts for quick and simple assembly.

## Easy to hold



The unique bottle shape makes this bottle easy to hold and grip in any direction.

## Easy to clean



With its wide neck and rounded corners, our bottle is easy to clean quickly and thoroughly. Ease of cleaning helps to ensure more hygienic feeding.

## Leak-free design

Our Anti-colic bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

## Growing up with your baby


Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

## BPA free



Philips Avent Anti-colic bottle is made of BPA free material (PP).

## Fits your baby's growing needs



The Philips Avent Anti-colic bottle range offers different teat flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All teats are available in twin packs: Newborn, Slow, Medium, Fast and Variable Flow, and Thick feed.

# Specifications

## Material

**Bottle:** BPA free, Polypropylene

**Teat:** BPA free, Silicone

## What is included

**Baby bottle:** 2 pcs

## Bottle

**Capacity:** 11oz/330ml

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Bottle use:** Dishwasher & microwave safe,  
Easy to assemble, Easy to clean, Easy to hold

## Development stages

**Stages:** 0-12 months

## Functions

**Anti-colic valve:** designed to reduce air ingestion

**Teat:** Easy latch on, Ribbed texture prevents teat collapse, Proven anti-colic system

**Ease of use:** easy to clean and assemble, leak-free design, 4 pieces for easy assembly



\* 0% BPA, following EU regulation 10/2011

\* At 2 weeks of age, babies fed with a Philips Avent bottle showed less colic, and significantly less fussing at night compared to babies fed with another competitor bottle.

\* Teat design proven to prevent teat collapse and associated air ingestion and feeding interruptions.

\* What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying and fussing.