

# PHILIPS

## AVENT

AVENT Combined  
Steamer and Blender



SCF870/21



## Steam, flip, blend and serve

### Feeding kids with Philips AVENT

You can easily prepare healthy homemade babyfood with the Philips AVENT Combined Steamer and Blender: first steam fruit, vegetables, fish or meat and then simply lift and flip the jar over to blend it, with no transfer of food required!

#### **Nutritious meals made easy**

- Compact design
- Easy water filling
- Steam and flip jar over to blend
- Few parts to clean

#### **Prepares babyfood for every stage of weaning**

- Blends also meat, fish and pulses
- Powerful blending performance
- Suitable for sauces and finger food

#### **Helps preserve nutrients & texture**

- Efficient steam cycle
- Retains cooking liquids

#### **Includes recipes & information on weaning**

- 12 age appropriate recipes
- Professional advice

# Highlights

**Compact design**

Compact design for easy storage.

**Blends also meat, fish and pul**

for complete meals

**Retains cooking liquids**

All cooking liquids are retained and can be reincorporated in blending

**12 age appropriate recipes**



for 3 stages of weaning

**Easy water filling**



Easy water filling

**Steam and flip jar over**



Steam and flip jar over to blend

**Suitable for sauces and finger**

Suitable for sauces and finger food

**Efficient steam cycle**



Efficient steam cycle

**Few parts to clean**

Few parts to clean

**Powerful blending performance**



for fine purees

**Professional advice**



from paediatric nutritionist and child psychologist

# Specifications

**Country of origin**

Turkey

**What is included**

- Steamer/blender: 1
- Spatula: 1
- Measuring cup: 1
- Recipe booklet: 1

**Technical specifications**

- Frequency: 50 - 60 Hz
- Power: 400 W
- Voltage: 220-240 V
- Safety Classification: Class 1
- Safety: Safety lock system for lid & bowl detection
- Capacity: 800 (solids) / 450 (liquids) ml
- Speed: 1
- Watertank volume: 200 ml
- Cord length: 70 cm
- Color/finishing: White/Green

**Weight and dimensions**

- F-box dimensions: 193 D x 243 W x 344 H mm
- Number of F-boxes in A-box: 2
- Product dimensions: 16.50 (round base) 30.8 (height) cm
- Product weight: 2 kg

**Material**

BPA free\*

**Development stages**

Stages: 6 months +, 1 year +, 6 - 12 months

