

PHILIPS

Screw ring for
feeding bottle

SCF925



Screw ring for bottle

Find matching products on the specifications tab

The adapter ring is one part of the Avent two-piece anti-colic system. This consists of the ring and a silicone teat. As your baby feeds, the unique skirt on the teat flexes to allow air into the feeding bottle instead of your baby's tummy.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types: SCF660/04, SCF660/07, SCF660/10, SCF660/17, SCF660/18, SCF660/20, SCF660/27, SCF660/28, SCF660/99, SCF663/00, SCF663/04, SCF663/07, SCF663/10, SCF663/17, SCF663/18, SCF663/20, SCF663/27, SCF663/28, SCF663/37, SCF663/47, SCF666/10, SCF666/17, SCF666/18, SCF666/27

