

PHILIPS

AVENT

Digital thermometer

Unique and playful design

Safe to play with

Waterproof and floats

Bath and bedroom measurement



SCH480/20



A smile for you and your baby

Happy moments from bath time to bed time

The Philips Avent digital thermometer is your bedroom and bath time buddy. It helps you easily monitor the room or water temperature while your little one inspects its playful petals.

Fast and reliable measurement

- Comfy bathing at a safe temperature
- Great sleep starts with the perfect room temperature
- Upright digital display

Safe to play with

- Tactile design with high safety standards
- Tactile design is safe and fun

Highlights

Bathe safely



Enjoy complete confidence in your baby's bath time comfort. Use the thermometer to check bath water stays between 36.5°C and 38°C. Keep it below 39°C for your baby's safety.

Sleep sweetly

Help your baby sleep soundly and sweetly by keeping the room temperature around 18°C. Not too cool. Not too cozy. Just right!

Upright digital display



We know you already have your hands full. So we designed our accurate digital display to sit upright for clear and easy reading, even at a glance.

Created for play

The bath and bedroom thermometer should be as fun for your baby as it is reassuring for you. So whatever the age of your baby, we made sure it complies with all the safety standards you would expect.

Tactile design is safe and fun

Fully waterproof and engineered to float, the bath and bedroom thermometer is a handy companion at bath time. So you and your baby can just enjoy another happy moment together.

Specifications

Technical specifications

Complies with: EN 71

Operating temperature range: 10–45 °C

Development stages

Stages: 0–12 months

