

Philips AVENT  
Digital baby thermometer  
set



SCH540/01

## Gentle and accurate

Combining the benefits of a digital thermometer and a soother thermometer.

### **Comfort for your child**

- Choice to adapt temperature measurement method to suit baby
- Soft, flexible tip digital thermometer
- Specially designed soother thermometer

### **Fast and reliable digital temperature readings**

- Clear measurement displays and battery status indicators
- Battery low indicator

### **Convenient for you and your baby**

- Includes AVENT Freeflow Soother

**PHILIPS**  
**AVENT**

# Specifications

## Country of origin

- Made in China: Thermometers, Yes
- Made in England: Soothers, Yes

## Technical specifications

- Measurement time: 10 secs  $\pm$ 3 secs in water (bath), 25 secs  $\pm$ 10 secs oral and rectal, 60 secs  $\pm$ 10 secs underarm
- Accuracy: 95–107.6  $\pm$ 0.2 °F
- Operating temperature range: 90–109.9 °F

## Accessories

- User manual

## Convenience

- Suitable for newborns
- Flexible tip
- Clear measurement display
- Automatic shut-off
- Battery status indicator

## Power

- Battery type: LR41
- Number of batteries: 1 per thermometer
- Operating time on battery: >100 hour(s)



Issue date 2019-07-30

Version: 4.0.1

12 NC: 0000 000 00000

© 2019 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)